



韦克斯科技（北京）有限公司

WIX TECHNOLOGY BEIJING CO., LTD

WIX-miniPRO4 Mini Vertical Electrophoresis Cell



1 ~ 4 piece(s) of manual gel(s) or premixed gel(s) is (are) available, 4 pieces of gel(s) electrophoresis can be made within 1 hour.

Features

- ◆ filler strip is fixed to glass plate to avoid leakage of gel.
- ◆ Gel-making kit ensures the accurate aligning.
- ◆ Excellent leakproofness.
- ◆ Special comb does not stop gel polymerization. Built-in kit isolates gel from air.
- ◆ Thickness and well-marked comb is readable.
- ◆ Equipped with multi-functional gel scraper.



韦克斯科技（北京）有限公司

WIX TECHNOLOGY BEIJING CO., LTD



Technical parameter

- ◆ Quantity: 1-4 piece(s)
- ◆ Thickness: 1.0mm: 10, 15 wells (standard), 0.75mm: 10 & 15 wells, 1.5mm: 10 & 15 wells (optional)
- ◆ Premixed gel: compatible with PIERCE, INVITROGEN, BIO-RAD
- ◆ Glass plate dimension: 100×83mm
- ◆ Gel dimension: 83×73mm
- ◆ Comb specification: 1.0mm thickness, 10 & 15 wells; 0.75mm, 10 & 15 wells; 1.5mm 10 & 15 wells.
- ◆ Dimension: 180×130×160mm (L×W×H)